

Softech solutions strongly adhere to compliance with MICROSOFT WINDOWS operating systems and specified hardware requirements herein. The entire Softech solution suite has been tested and certified using the system requirements described herein. Any deviation may result in unpredictable behavior and/or poor performance for which Softech is not liable.

<p><b>Windows Server for Dentech and DenChart</b></p>	<p><b><u>Hardware</u></b></p> <ul style="list-style-type: none"> <li>• Core 2 Duo, i5, Xeon Dual Core, Xeon Quad Core 32-bit processor</li> <li>• 4 GB RAM</li> <li>• 500 GB hard drive</li> <li>• DVD-ROM drive</li> <li>• Ethernet 1GB network</li> <li>• LCD Monitor</li> <li>• Mouse/keyboard</li> <li>• 2 USB 2.0 capable ports</li> </ul> <p>UPS/Surge suppression hardware</p> <p>Backup should be adequate for backing up the entire server hard disk including the operating system (full system backup)</p> <p>High speed Internet access is highly recommended (3mb/sec DL - 1mb/sec UP speed minimum)</p> <p>Trend Virus protection with updates is highly recommended with exclusions enabled</p>	<p><b><u>Supported Operating Systems</u></b></p> <p><b><u>Single User</u></b> Dentech software multi-user version Windows XP Professional SP3 32-bit Windows Vista Business or higher 32-bit** Windows 7 Professional 32-bit</p> <p><b><u>With 2-4 Workstations</u></b> Dentech software multi-user version Peer-to-Peer Configuration TCP/IP Protocol Windows XP Professional SP3 32-bit Windows Vista Business or higher 32-bit** Windows 7 Professional 32-bit</p> <p><b><u>With 5+ Workstations</u></b> Dentech software multi-user version Server/Workstation Configuration using dedicated server TCP/IP Protocol Windows Server 2003 or 2008 32-bit</p>
<p><b>Windows Workstation for Clinical &amp; Business Operations (including digital imaging)</b></p>	<p><b><u>Hardware</u></b></p> <ul style="list-style-type: none"> <li>• Pentium Dual Core, Core 2 Duo, i3, i5, i7 32 bit processor</li> <li>• 3 GB RAM(XP); 4 GB (Win7)</li> <li>• 80 GB available disk space (New equipment –160GB recommended)</li> <li>• DVD-ROM drive</li> <li>• Ethernet 1GB network</li> <li>• 1024x768, 32 bit capable video card with 256 MB RAM or higher</li> <li>• Standard LCD monitor with a minimum of 1024x768 screen resolution (Dual monitors are recommended)</li> <li>• USB Mouse/keyboard</li> <li>• 2 or more powered USB 2.0 ports (for imaging 4 ports are recommended) 1 Dedicated port for headset (For voice enabled products)</li> <li>• Additional PCI Express expansion slots may be required</li> <li>• DenChart certified headset*</li> </ul> <p>UPS/Surge suppression hardware</p> <p>High speed Internet access is highly recommended (3mb/sec DL - 1mb/sec UP speed minimum) Trend Virus protection with updates is highly recommended with exclusions enabled</p>	<p><b><u>Supported Operating Systems</u></b></p> <p>Windows XP Professional SP3 32-bit Windows 7 Professional 32-bit (with Dentech 8.0 and higher installed)</p> <p><i>With -</i> Crystal Reports 2008 Windows Installer 3.1 Windows .Net 3.5 SAPI 5.1 Engine (XP)/8.0 Engine (Win7) * Wingdings II Font</p>

\* Not required on Business Workstations.

\*\* End of Life January 1, 2012