

## 1 Specifications

Model	1001-IOM	
Camera	Power Source	3.7V
	Size	Approx. 266 x 31 x 37 mm
	Weight	Approx. 140g
	Resolution	1.3 Mega pixel sensor
	Sensor/Lens Type	CMOS / Liquid Lens
Charger	Measurement	Approx. 169 x 148 x 89 mm
	Weight	Approx. 170g
	Operation Voltage	5V
Disposable cover	Dimensions	Approx. 32 x 100 mm
	Quantity	Approx. 100 pcs by weight x 2 ea
	Raw material	Polyethylene
Adapter	Input Voltage	AC 100 ~ 240 V
	Output Voltage	DC 5V 2A
	UL Number	E333054

## 2 Minimum PC System

PC system:	2.0 GHz CPU, 1GB RAM
OS:	Windows XP service pack 2 (Windows XP 64 bit is not supported)
USB Port:	USB 2.0



## 1 Specifications

Model	1001-IOM	
Camera	Power Source	5V (USB Power)
	Size	Approx. 202 x 30 x 28 mm
	Weight	Approx. 56g
	Resolution	1.3 Mega pixel sensor
	Sensor/Lens Type	CMOS / Liquid Lens
USB cable	Cable Length	Approx. 3 M (9.8 ft)
Disposable cover	Quantity	Approx. 100 pcs by weight x 2 ea
	Raw material	Polyethylene

## 2 Minimum PC System

PC system:	2.0 GHz CPU, 1GB RAM
OS:	Windows XP service pack 2
USB Port:	USB 2.0

